***** usage instruction for MPMS since 2007.08.01 *****

before your measurement

1. fill up the log note below. *¹
2. run "log MPMS". the file name must be noted in the log note below.
(*Utilities/Log MPMS Data/*, fill and check the required parameters, log-file name, sampling rate (~100 sec), temperature, field etc.. in the dialog box, and push the "START" button)
3. write any anomalous thing you feel.

in your measurement

1. keep the He-level of MPMS above 50% when you set $H > 1$ T.
2. check the liq-N₂ in the outer jacket of MPMS. if there is no liq-N₂, please refill it. (be sure the exhaust valve sequence when you do liq-N₂ filling)
3. keep the amount of liq-He in the vessel more than 10 liter.
4. return the vessel to the cryogenics service laboratory and/or contact to me as soon as possible when the residual liq-He is less than 10 liter.

after your measurement

1. set $H = 0$ Oe.
2. finish the "log MPMS".
3. turn the MPMS standby mode for economical consuming of liq-He.
(*Instrument/Standby*, carefully see the dialog box and push "OK" only one time)
4. convert MPMS to the default setting (DC-setting).
5. do not leave your samples.
6. clean up the desk and the around it to make these more clear than that before you use.

Touru Yamauchi
emag@issp.u-tokyo.ac.jp
04-7136-3436
090-1201-3643 (emergency only)

***** MPMS log note *****

individual	
* name	
* affiliation	
* phone number (mobile phone is favorable)	
* usage term	yy/mm/dd : / / ~ / /
* sample	
* file name of "log MPMS"	
measurements	
* option (setting)	DC / RSO / High Temp
* probe	DC-1 / DC-2 / RSO-1 / RSO-2
* probe status	injured / uninjured
* temperature & field range	$\leq T \leq$ (K), $\leq H \leq$ (T)
He fill ; if you did	
vessel number	
residual amount of liq-He	(liter)
after filling, He level in MPMS	(%)
all about your measurements	
remarks	

¹ asterisk (*) must be filled. only one blank makes you difficult to get your machine time of MPMS

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